

PURESITE - Closed Male Connector


B | BRAUN

GB	Instructions for use	7
DE	Gebrauchsanweisung	8
BG	Инструкция за употреба	9
CN	使用说明	10
CZ	Návod k použití	11
DK	Bruksanvisning	12
EE	Kasutusjuhend	13
ES	Instrucciones de uso	14
FI	Käyttöohje	15
FR	Mode d'emploi	16
GR	Οδηγίες χρήσης	17
HR	Upute za korištenje	18
HU	Használati útastítás	19
IT	Istruzioni per l'uso	20
LT	Naudojimo instrukcija	21
LV	Lietošanas instrukcija	22
NL	Gebruiksaanwijzing	23
NO	Bruksanvisning	24
PL	Instrukcja użytkowania	25
PT	Instruções de utilização	26
RO	Instrucțiuni de utilizare	27
RU	Инструкция по применению	28
SE	Bruksanvisning	29
SI	Navodila za uporabo	30
SK	Návod na použitie	31
TR	Kullanım Kılavuzu	32


I. A. Handling PureSite - Closed Male Connector

I. B. + II. B. Handling PureSite - Closed Male Connector

II. A. Handling PureSite - Closed Male Connector

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Please read carefully!

Materials used: PSU (Polysulfone) ; PC (Polycarbonate); Silicone

CAUTION: Ensure medications are prepared and administered in compliance with drug manufacturer's recommendations.

- PureSite is designed for use with luer lock syringes and access devices.
- Single use only, Do not resterilize. Use only if the packaging is intact.
- PureSite is not compatible with the solvent N,N-Dimethylacetamide commonly found in Busulfex®/Busulfan.
- Disposal of components should be in accordance with local regulations and procedures for the handling of hazardous drugs
- Use Aseptic technique (only alcohol pads).
- Swabbing as per manufacturer and institution's protocol

I. A. For Medication Preparation (Admixture) using a syringe:

1. Open package of PureSite at designated end. Attach female luer end of PureSite to syringe by holding the orange grip with the package (1a).
2. The female luer of PureSite has a permanent connection and cannot be detached from the syringe (spinning) (1b).
3. Connect PureSite to needle free vial access device (2).
4. Invert vial and aspirate medication.
5. To expel air from syringe, draw half of required volume, tap inverted syringe. Depress plunger to expel air, then draw up plunger to total required volume (3).
6. Disconnect PureSite from needle free vial access device.
7. Cap for better protection during transportation (See dust cover cap's instructions) (4).
8. Transport per institutional protocol.

I. B. For IV Bolus Administration using Syringe:

1. Swab luer access device on Patient IV Administration set per institutional protocol.

2. Swab PureSite male luer (5) and firmly connect the PureSite with luer access device. Hold the PureSite (orange grip) instead of the syringe (6).
3. Deliver medication.
4. Disconnect by removing PureSite (orange grip) from luer access device.
5. Dispose of syringe assembly per institutional protocol.

II. A. For Medication Preparation using a Secondary IV Administration set:

1. Prepare IV solution container and secondary administration set per manufacturer's instructions.
2. Attach secondary set to solution container and prime tubing per manufacturer's instructions.
3. Open package of PureSite at designated end. Attach female luer end of PureSite to end of Secondary IV Administration set by holding the orange grip with the package (7a).
4. The female luer of PureSite has a permanent connection and cannot be detached from the IV set (spinning) (7b).
5. Remove residual air from PureSite by connecting a female luer connector (opening of PureSite) and priming (8).
6. Remove female luer connector.
7. Cap if desired (See dust cover cap's instructions) (9).
8. Transport per institutional protocol.

II. B. For Medication IV Administration using a Secondary IV Administration set:

1. Lower patient's primary IV container and hang secondary IV administration set per manufacturer's instructions.
2. Swab luer access device on patient IV Administration set.
3. Remove cap (if used) from male luer of PureSite. Swab PureSite (5) and attach to female luer of access device (hold the PureSite with the orange grip) (6).
4. Proceed with fluid administration
5. Remove secondary IV set by disconnecting PureSite from luer access device on primary IV line.
6. Dispose of secondary IV set per institutional protocol.

Warning

Re-use of single-use devices creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability. Contamination and/or limited functionality of the device may lead to injury, illness or death of the patient.

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B. Braun Melsungen AG